


Progesterone Receptor Mouse mAb

ABZ4064


Product Overview

Name	Progesterone Receptor Mouse mAb
Accession(Primary)	P06401
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat
Recommended Dilution	IHC-1:50-1:100; IF-1:50-1:200


Target information (P06401)

Gene Name PGR**About Gene Symbol**


Immunohistochemistry analysis of paraffin-embedded Human uterus tissue using Progesterone Receptor (8H3) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Negative control was used by secondary antibody only.




Immunohistochemical analysis of paraffin-embedded Human tonsils using Progesterone Receptor (8H3) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Negative control was used by secondary antibody only.




Immunohistochemistry analysis of paraffin-embedded mouse testis tissue using PR antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Negative control was used by secondary antibody only.



Immunofluorescence analysis of Progesterone Receptor (8H3) in Human appendix tissue using Progesterone Receptor (8H3) antibody (Z15) (red), and DAPI (blue).



Immunofluorescence analysis of Progesterone Receptor (8H3) in mouse liver using PR antibody (Z15) (red), and DAPI (blue).



Immunofluorescence analysis of Progesterone Receptor in rat heart using Progesterone Receptor (8H3) antibody (red) and DAPI (blue).