

CD1a (4C3) Mouse mAb

ABZ4065


Product Overview

Name	CD1a (4C3) Mouse mAb
Accession(Primary)	P06126
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat
Recommended Dilution	IHC-1:50-1:100 ; IF-1:50-1:200


Target information (P06126)

Gene Name CD1A


About Gene Symbol




Immunohistochemistry analysis of paraffin-embedded Human Tonsil tissue using CD1a (4C3) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Negative control was used by secondary antibody only.




Immunohistochemical analysis of paraffin-embedded Human tonsils using CD1a (4C3) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Negative control was used by secondary antibody only.



Immunohistochemistry analysis of paraffin-embedded mouse heart tissue using CD1a antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Negative control was used by secondary antibody only.



Immunofluorescence analysis of CD1a (4C3) in mouse heart tissue using CD1a (4C3) antibody (9H6) (red), and DAPI (blue).



Immunohistochemistry analysis of paraffin-embedded Human tonsils using CD1a (4C3) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.