

CD2 (3F10) Mouse mAb

ABZ4066


Product Overview

Name	CD2 (3F10) Mouse mAb
Accession(Primary)	P06729
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat
Recommended Dilution	IHC-1:50-1:100


Target information (P06729)

Gene Name CD2


About Gene Symbol




Immunohistochemistry analysis of paraffin-embedded Human colon tissue using CD2 (3F10) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Negative control was used by secondary antibody only.



Immunohistochemical analysis of paraffin-embedded Human tonsils using CD2 (3F10) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Negative control was used by secondary antibody only.



Immunohistochemistry analysis of paraffin-embedded mouse spleen tissue using CD2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Negative control was used by secondary antibody only.



Immunohistochemistry analysis of paraffin-embedded Human Tonsil Carcinoma using CD2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.