


SGK1 Mouse mAb

ABZ4071


Product Overview

Name	SGK1 Mouse mAb
Accession(Primary)	O00141
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat
Recommended Dilution	IHC-1:50-1:100

Target information (O00141)

Gene Name SGK1**About Gene Symbol**

Immunohistochemistry analysis of paraffin-embedded Human Brain Tissue using SGK1 (6E4) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemical analysis of paraffin-embedded Human tonsils using SGK1 (6E4) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.