

## Cystatin C (7E3) Mouse mAb


ABZ4073

### Product Overview


|                             |                                                |
|-----------------------------|------------------------------------------------|
| <b>Name</b>                 | Cystatin C (7E3) Mouse mAb                     |
| <b>Accession(Primary)</b>   | P01034                                         |
| <b>Host Species</b>         | Mouse                                          |
| <b>Clonality</b>            | Monoclonal                                     |
| <b>Isotype</b>              | IgG1                                           |
| <b>Species Reactivity</b>   | Human                                          |
| <b>Calculated MW</b>        | Calculated MW: 16 kDa; Observed MW: 16 kDa     |
| <b>Recommended Dilution</b> | WB-1:500-1:1000; IHC-1:50-1:100; ELISA-1:10000 |

### Target information (P01034)


|                          |                                                                                                                                                                                                                                                                                                    |
|--------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | CST3                                                                                                                                                                                                                                                                                               |
| <b>About Gene Symbol</b> | CST3 (Cystatin C) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets. |



Immunohistochemistry analysis of paraffin-embedded Human Kidney Tissue using Cystatin C (7E3) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemical analysis of paraffin-embedded Human tonsils using Cystatin C (7E3) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human Brain Tissue using Cystatin C antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.