

PDGFR alpha (7H6) Mouse mAb

ABZ4080


Product Overview

Name	PDGFR alpha (7H6) Mouse mAb
Accession(Primary)	P16234
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat
Recommended Dilution	IHC-1:50-1:100

Target information (P16234)

Gene Name PDGFRA

About Gene Symbol



Immunohistochemistry analysis of paraffin-embedded rat Skin Tissue using PDGFR alpha (7H6) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemical analysis of paraffin-embedded Human tonsils using PDGFR alpha (7H6) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.