

c-Jun Mouse mAb


ABZ4081

Product Overview


Name	c-Jun Mouse mAb
Accession(Primary)	P05412
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat
Recommended Dilution	IHC-1:50-1:100

Target information (P05412)

Gene Name	JUN
About Gene Symbol	JUN (Jun Proto-Oncogene, AP-1 Transcription Factor Subunit) is a proto-oncogene involved in gene transcription and regulation. It is an important research target in cancer research and infectious disease and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Immunohistochemistry analysis of paraffin-embedded Human Colon Carcinoma Tissue using c-Jun (10G2) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemical analysis of paraffin-embedded Human tonsils using c-Jun (10G2) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.