

TATA binding protein TBP Mouse mAb

ABZ4084


Product Overview

Name	TATA binding protein TBP Mouse mAb
Accession(Primary)	P20226
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat
Recommended Dilution	IHC-1:50-1:100

Target information (P20226)

Gene Name TBP

About Gene Symbol



Immunohistochemistry analysis of paraffin-embedded Human Breast Carcinoma using TATA Box Binding Protein (10C1) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.