

S6K1 (3A2) Mouse mAb

ABZ4085


Product Overview

Name	S6K1 (3A2) Mouse mAb
Accession(Primary)	P23443
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat
Recommended Dilution	IHC-1:50-1:100


Target information (P23443)

Gene Name RPS6KB1

About Gene Symbol



Immunohistochemistry analysis of paraffin-embedded Human Lung Carcinoma Tissue using S6K1 (3A2) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemical analysis of paraffin-embedded Human tonsils using S6K1 (3A2) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.