

## JAK1 (8D8) Mouse mAb


ABZ4086

### Product Overview


|                             |                      |
|-----------------------------|----------------------|
| <b>Name</b>                 | JAK1 (8D8) Mouse mAb |
| <b>Accession(Primary)</b>   | P23458               |
| <b>Host Species</b>         | Mouse                |
| <b>Clonality</b>            | Monoclonal           |
| <b>Isotype</b>              | IgG1                 |
| <b>Species Reactivity</b>   | Human, Mouse, Rat    |
| <b>Recommended Dilution</b> | IHC-1:50-1:100       |

### Target information (P23458)

|                          |                                                                                                                                                                                                                                                                                          |
|--------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | JAK1                                                                                                                                                                                                                                                                                     |
| <b>About Gene Symbol</b> | JAK1 (Janus Kinase 1) is a receptor tyrosine kinase involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets. |



Immunohistochemistry analysis of paraffin-embedded Human Lung Carcinoma Tissue using JAK1 (8D8) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemical analysis of paraffin-embedded Human tonsils using JAK1 (8D8) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.