

PHC1 (8G5) Mouse mAb

ABZ4088


Product Overview

Name	PHC1 (8G5) Mouse mAb
Accession(Primary)	P78364
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human
Recommended Dilution	WB-1:500-1:1000; FC-1:100


Target information (P78364)

Gene Name PHC1

About Gene Symbol



Western blot analysis of HEK293 (1) and PHC1 (AA: 7581004)hIgGFc transfected HEK293 (2) lysates using PHC1 antibody.



Flow cytometry analysis of HEK293 stained with PHC1 antibody (green) and negative control (red).