

SMN/Gemin 1 (2F1) Mouse mAb


ABZ4089

Product Overview


Name	SMN/Gemin 1 (2F1) Mouse mAb
Accession(Primary)	Q16637
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Cat
Calculated MW	Calculated MW: 32 kDa; Observed MW: 39 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IP-1:20; FC-1:50-1:100

Target information (Q16637)

Gene Name	SMN1
About Gene Symbol	SMN1 (Survival Of Motor Neuron 1, Telomeric) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Gemin 1 in COS7, Jurkat and HeLa lysates using Gemin 1 antibody.



Immunoprecipitation analysis of SMN (2F1) in HeLa lysates using Gemin 1 antibody.