

Prohibitin (5H7) Mouse mAb


ABZ4090

Product Overview

Name	Prohibitin (5H7) Mouse mAb
Accession(Primary)	P35232
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat, Cat
Calculated MW	Calculated MW: 30 kDa; Observed MW: 30 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200; FC-1:50-1:100

Target information (P35232)

Gene Name	PHB1
About Gene Symbol	PHB1 (Prohibitin 1) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Prohibitin in C6, COS7, 3T3 and HeLa lysates using Prohibitin antibody.