

Dynamin 1 (5B9) Mouse mAb


ABZ4091

Product Overview


Name	Dynamin 1 (5B9) Mouse mAb
Accession(Primary)	Q05193
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Species Reactivity	Human
Calculated MW	Calculated MW: 97 kDa; Observed MW: 97 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100

Target information (Q05193)


Gene Name	DNM1
About Gene Symbol	DNM1 (Dynamin 1) is a dynamin GTPase involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Western blot analysis of Dynamin 1 (5B9) in C6 , NIH/3T3 , SK-N-SH , LN18 , SH-SY5Y lysates and rat brain lysates using Dynamin 1 (5B9) antibody.



Immunohistochemical analysis of paraffin-embedded Human tonsils using Dynamin 1 (5B9) antibody, with DAB staining. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



B



C