

RON (9B6) Mouse mAb

ABZ4093


Product Overview

Name	RON (9B6) Mouse mAb
Accession(Primary)	Q04912
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human
Calculated MW	Calculated MW: 152 kDa; Observed MW: 152 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (Q04912)

Gene Name MST1R

About Gene Symbol



Western blot analysis of RON (9B6) in HCC827, HT-29, HCT116 and BxPC-3 lysates using RON (9B6) antibody.

