

ApoA4 Mouse mAb

ABZ4094


Product Overview

Name	ApoA4 Mouse mAb
Accession(Primary)	P06727
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Species Reactivity	Human
Calculated MW	Calculated MW: 45 kDa; Observed MW: 45 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P06727)

Gene Name APOA4

About Gene Symbol



Western blot analysis of ApoA4 in Human serum , Human plasma , HepG2 lysates and SMMC7721 lysates using ApoA4 antibody.

