

Prostate Specific Antigen Mouse mAb

ABZ4095


Product Overview

Name	Prostate Specific Antigen Mouse mAb
Accession(Primary)	P07288
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human
Calculated MW	Calculated MW: 29 kDa; Observed MW: 37 kDa
Recommended Dilution	WB-1:500-1:1000; FC-1:50-1:100


Target information (P07288)

Gene Name KLK3**About Gene Symbol**

Human KLK3 Recombinant Protein



Western blot analysis of Prostate Specific Antigen in Human KLK3 (AA: 10921recombinant protein.) using Prostate Specific Antigen antibody.



Flow Cytometry analysis of HeLa cells stained with KLK3 antibody (green) and negative control (red).