

Cytokeratin 8 (10A5) Mouse mAb

ABZ4096


Product Overview

Name	Cytokeratin 8 (10A5) Mouse mAb
Accession(Primary)	P05787
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 54 kDa; Observed MW: 54 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200


Target information (P05787)

Gene Name KRT8


About Gene Symbol




Immunohistochemistry analysis of paraffin-embedded Human colon tissue using Cytokeratin 8 (10A5) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Negative control was used by secondary antibody only.




Immunohistochemical analysis of paraffin-embedded Human tonsils using Cytokeratin 8 (10A5) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Negative control was used by secondary antibody only.



Immunohistochemistry analysis of paraffin-embedded mouse liver tissue using CK8 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Negative control was used by secondary antibody only.



Western blot analysis of Cytokeratin 8 (10A5) in HeLa, mouse Brain lysates using Cytokeratin 8 (10A5) antibody.



Immunohistochemistry analysis of paraffin-embedded human colon cancer tissue using Cytokeratin 8 (10A5) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.