

Cytokeratin 6 (1A5) Mouse mAb

ABZ4097


Product Overview

Name	Cytokeratin 6 (1A5) Mouse mAb
Accession(Primary)	P02538
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human
Calculated MW	Calculated MW: 60 kDa; Observed MW: 56-60 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200

Target information (P02538)

Gene Name KRT6A/KRT6B/KRT6C

About Gene Symbol



Immunohistochemistry analysis of paraffin-embedded Human breastcancer tissue using Cytokeratin 6 (1A5) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Negative control was used by secondary antibody only.

