


CD16 Mouse mAb

ABZ4099


Product Overview

Name	CD16 Mouse mAb
Accession(Primary)	P08637
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human
Calculated MW	Calculated MW: 29 kDa; Observed MW: 45 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100

Target information (P08637)

Gene Name FCGR3A**About Gene Symbol**

Immunohistochemistry analysis of paraffin-embedded Human liver tissue using CD16 (9G6) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Negative control was used by secondary antibody only.



Western blot analysis of CD16 (9G6) in Jurkat, K562 lysates using CD16 (9G6) antibody