

Hsp90 beta (6A7) Mouse mAb


ABZ4100

Product Overview

Name	Hsp90 beta (6A7) Mouse mAb
Accession(Primary)	P08238
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 83 kDa; Observed MW: 90 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200

Target information (P08238)

Gene Name	HSP90AB1
About Gene Symbol	HSP90AB1 (Heat Shock Protein 90 Alpha Family Class B Member 1) is a molecular chaperone involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




A

B

C

Immunofluorescence analysis of Hsp90 beta (6A7) in mouse spleen tissue using Hsp90 beta (6A7) antibody(red),and DAPI (blue).




A

B

C

Immunofluorescence analysis of Hsp90 beta (6A7) in rat lung tissue using HSP90 β antibody(red),and DAPI (blue).



Western blot analysis of Hsp90 beta (6A7) in HeLa, mouse Brain tissue, rat Brain tissue lysates using Hsp90 beta (6A7) antibody