


Rabbit IgG H&L (10D2) Mouse mAb

ABZ4101

Product Overview

Name	Rabbit IgG H&L (10D2) Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Rabbit
Calculated MW	Calculated MW: 25,55 kDa; Observed MW: 25,55 kDa
Recommended Dilution	WB-1:500-1:1000



Western blot analysis of Reduced Rabbit IgG, H+L using Rabbit IgG H&L (10D2) antibody.