


Rabbit IgG Fc (6D3) Mouse mAb

ABZ4102

Product Overview

Name	Rabbit IgG Fc (6D3) Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Rabbit
Calculated MW	Calculated MW: 55 kDa; Observed MW: 55 kDa
Recommended Dilution	WB-1:500-1:1000



Western blot analysis of Reduced Rabbit IgG, H+L using Rabbit IgG Fc (6D3) antibody.