


Rabbit IgG Light Chain (7F9) Mouse mAb

ABZ4103

Product Overview

Name	Rabbit IgG Light Chain (7F9) Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Rabbit
Calculated MW	Calculated MW: 25 kDa; Observed MW: 25 kDa
Recommended Dilution	WB-1:500-1:1000



Western blot analysis of Reduced Rabbit IgG, H+L using Rabbit IgG Light Chain (7F9) antibody.