

Ubiquitin (1B8) Mouse mAb

ABZ4106


Product Overview

Name	Ubiquitin (1B8) Mouse mAb
Accession(Primary)	P0CG47
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Refer to figures
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100


Target information (P0CG47)

Gene Name UBB


About Gene Symbol



Western blot analysis of Ubiquitin (1B8) in HeLa lysates, 3T3 lysates, rat Brain lysates using Ubiquitin antibody.



Immunohistochemistry analysis of paraffin-embedded Human Breast Carcinoma Tissue using Ubiquitin antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human Stomach Carcinoma Tissue using Ubiquitin antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.