


HP1 gamma (9E5) Mouse mAb

ABZ4109


Product Overview

Name	HP1 gamma (9E5) Mouse mAb
Accession(Primary)	Q13185
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 21 kDa; Observed MW: 24 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100


Target information (Q13185)

Gene Name CBX3**About Gene Symbol**


Immunohistochemistry analysis of paraffin-embedded Human Colon Carcinoma Tissue using HP1 gamma (9E5) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemical analysis of paraffin-embedded Human tonsils using HP1 gamma (9E5) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human Colon Carcinoma Tissue using HP1 gamma (9E5) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of HP1 gamma (9E5) in HeLa lysates, 3T3 lysates, PC-12 lysates using HP1 gamma (9E5) antibody.