

Hsc70 (5C2) Mouse mAb


ABZ4112

Product Overview


Name	Hsc70 (5C2) Mouse mAb
Accession(Primary)	P11142
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 71 kDa; Observed MW: 71 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100

Target information (P11142)

Gene Name	HSPA8
About Gene Symbol	HSPA8 (Heat Shock Protein Family A (Hsp70) Member 8) is a molecular chaperone involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mecha



Immunohistochemistry analysis of paraffin-embedded Human Breast Carcinoma using Hsc70 (5C2) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of Hsc70 (5C2) in HeLa, A549, HL-60, U2OS, C6 lysates using Hsc70 (5C2) antibody