


TATA binding protein TBP Mouse mAb

ABZ4114

Product Overview

Name	TATA binding protein TBP Mouse mAb
Accession(Primary)	P20226
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 38 kDa; Observed MW: 35-45 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P20226)

Gene Name TBP**About Gene Symbol**

Western blot analysis of TATA Box Binding Protein (9D1) in HeLa, mouse Brain, rat Brain lysates using TBP/TATA Binding Protein antibody.

