

TNF alpha (5F3) Mouse mAb

ABZ4115


Product Overview

Name	TNF alpha (5F3) Mouse mAb
Accession(Primary)	P01375
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human
Recommended Dilution	ELISA-1:10000

Target information (P01375)

Gene Name

TNF

About Gene Symbol

Standard Curve for TNFa: Capture Antibody mouse mAb 200946 to TNFa at 4ug/ml and AntiTNF Human antibody for detecting.