


Ferritin Light Chain (1A2) Mouse mAb

ABZ4116

Product Overview

Name	Ferritin Light Chain (1A2) Mouse mAb
Accession(Primary)	P02792
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Species Reactivity	Human
Recommended Dilution	ELISA-1:10000

Target information (P02792)

Gene Name FTL**About Gene Symbol**

Indirect ELISA assay for antiFerritin Light Chain mouse mAb. Antigen coating concentration: 4ug/ml.