


Human IgG (2D7) Mouse mAb


ABZ4117

Product Overview

Name	Human IgG (2D7) Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human
Recommended Dilution	ELISA-1:10000



Recombinant Anthrax Protective Antigen used as Coating antigen, followed by AntiAnthrax Protective Antigen Human antibody and then incubated with AntiHuman IgG(201068) mouse mAb, HRP conjugated goat anti mouse IgG was used for coloring.



Recombinant Anthrax Protective Antigen used as Coating antigen, followed by different dilutions of supernatants which includes secreted AntiAnthrax Protective Antigen Human antibody, HRP conjugated AntiHuman IgG(201068) mouse mAb was used for coloring.