


Anthrax Protective gen (5C1) Mouse mAb

ABZ4118

Product Overview

Name	Anthrax Protective gen (5C1) Mouse mAb
Accession(Primary)	P13423
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Species Reactivity	Bacillus anthracis
Recommended Dilution	ELISA-1:10000

Target information (P13423)

Gene Name pagA**About Gene Symbol**

Standard Curve for Anthrax Protective Antigen:
Capture Antibody mouse mAb 201240 to Anthrax
Protective Antigen at 4ug/ml and other antibody for
detecting.