

FABP3 (3E8B10) Mouse mAb

ABZ4122


Product Overview

Name	FABP3 (3E8B10) Mouse mAb
Accession(Primary)	P05413
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Species Reactivity	Human
Recommended Dilution	ELISA-1:10000

Target information (P05413)

Gene Name FABP3

About Gene Symbol



Standard Curve for HFABP: Capture Antibody mouse mAb (5C4F5) to HFABP at 4ug/ml and Detector Antibody mouse mAb(3E8B10u30017B10G11u30014H1B8)to HFABP at 0.08 ug/ml.