

beta 2 Microglobulin (4C3) Mouse mAb


ABZ4124

Product Overview

Name	beta 2 Microglobulin (4C3) Mouse mAb
Accession(Primary)	P61769
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Species Reactivity	Human
Recommended Dilution	ELISA-1:10000

Target information (P61769)

Gene Name	B2M
About Gene Symbol	B2M (Beta-2-Microglobulin) is a membrane receptor involved in adaptive immunity. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



sELISA standard Curve for u03b22microglobulin.