

Lipocalin 2 (3C6) Mouse mAb

ABZ4129


Product Overview

Name	Lipocalin 2 (3C6) Mouse mAb
Accession(Primary)	P80188
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Species Reactivity	Human, Rat
Recommended Dilution	ELISA-1:10000

Target information (P80188)

Gene Name LCN2

About Gene Symbol



AntiFITC mouse mAb was used to capture FITClabeled NGAL mouse mAb (2011793C6), HRPlabeled NGAL mouse mAb (2011792F6) was used for detecting.