


Lipocalin 2 (4C12) Mouse mAb

ABZ4130

Product Overview

Name	Lipocalin 2 (4C12) Mouse mAb
Accession(Primary)	P80188
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Species Reactivity	Human
Recommended Dilution	ELISA-1:10000

Target information (P80188)

Gene Name LCN2**About Gene Symbol**

Indirect ELISA assay for antiNGAL mouse mAb. Antigen coating concentration: 4ug/ml.