

Creatine Kinase MM (6A12) Mouse mAb

ABZ4133


Product Overview

Name	Creatine Kinase MM (6A12) Mouse mAb
Accession(Primary)	P06732
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Species Reactivity	Human
Recommended Dilution	ELISA-1:10000

Target information (P06732)

Gene Name CKM

About Gene Symbol



Indirect ELISA assay for antiCKMM mouse mAb. CKBB as control antigen. Antigen coating concentration: 4ug/ml.