

BNP(3C9) Mouse mAb

ABZ4139


Product Overview

Name	BNP(3C9) Mouse mAb
Accession(Primary)	P16860
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human
Recommended Dilution	ELISA-1:10000


Target information (P16860)

Gene Name NPPB

About Gene Symbol



sELISA standard Curve for NTproBNP: Capture Antibody mouse mAb (2012893C9) to NTproBNP and Detector Antibody mouse mAb(2012891B5) to NTproBNP.



AntiFITC mouse mAb was used to capture
FITClabeled NTproBNP mouse mAb (2012893C9),
HRPlabeled Natriuretic Peptides B mouse mAb
(2012891B5) was used for detecting.