

BNP(5G5) Mouse mAb

ABZ4140


Product Overview

Name	BNP(5G5) Mouse mAb
Accession(Primary)	P16860
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human
Recommended Dilution	ELISA-1:10000


Target information (P16860)

Gene Name NPPB

About Gene Symbol



Indirect ELISA assay for antiNTproBNP mouse mAb. Antigen coating concentration: 4ug/ml.

[illegible]