


GFAP (9A2) Mouse mAb

ABZ4143

Product Overview

Name	GFAP (9A2) Mouse mAb
Accession(Primary)	P14136
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 50 kDa; Observed MW: 50 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200

Target information (P14136)

Gene Name GFAP**About Gene Symbol**

Immunohistochemistry analysis of paraffin-embedded Human liver tissue using GFAP (9A2) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Negative control was used by secondary antibody only.

