

CEACAM5 (2A8) Mouse mAb

ABZ4144


Product Overview

Name	CEACAM5 (2A8) Mouse mAb
Accession(Primary)	P06731
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human
Recommended Dilution	IHC-1:50-1:100; IF-1:50-1:200


Target information (P06731)

Gene Name CEACAM5


About Gene Symbol



Immunohistochemistry analysis of paraffin-embedded Human liver tissue using Carcino Embryonic Antigen CEA (2A8) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Negative control was used by secondary antibody only.



Immunofluorescence analysis of Carcino Embryonic Antigen CEA (2A8) in Human lung cancer tissue using Carcino Embryonic Antigen CEA (2A8) antibody (red), and DAPI (blue).



Immunohistochemistry analysis of paraffin-embedded Human colon cancer tissue using Carcino Embryonic Antigen CEA (2A8) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.