


TBL1X (4H2) Mouse mAb

ABZ4146


Product Overview

Name	TBL1X (4H2) Mouse mAb
Accession(Primary)	O60907
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Species Reactivity	Human
Calculated MW	Calculated MW: 62 kDa; Observed MW: 58 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200

Target information (O60907)

Gene Name TBL1X**About Gene Symbol**

Western blot analysis of TBL1X in SW480 lysates using TBL1X antibody.



Immunocytochemistry analysis of TBL1X (4H2) in HeLa using TBL1X antibody.