


Nucleophosmin Mouse mAb

ABZ4147

Product Overview

Name	Nucleophosmin Mouse mAb
Accession(Primary)	P06748
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 33 kDa; Observed MW: 38 kDa
Recommended Dilution	WB-1:500-1:1000; ChIP-1:20

Target information (P06748)

Gene Name NPM1**About Gene Symbol**

Western blot analysis of NPM1 in HeLa, Jurkat and 3T3 lysates using NPM1 antibody.

