

KLF4 Mouse mAb

ABZ4148

Product Overview

Name	KLF4 Mouse mAb
Accession(Primary)	O43474
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 55 kDa; Observed MW: 55,65 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (O43474)

Gene Name	KLF4
About Gene Symbol	KLF4 (KLF Transcription Factor 4) is a transcription factor involved in angiogenesis. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

