


Oct-4 (3F5) Mouse mAb

ABZ4149


Product Overview

Name	Oct-4 (3F5) Mouse mAb
Accession(Primary)	P20263
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Species Reactivity	Mouse
Calculated MW	Calculated MW: 38 kDa; Observed MW: 45 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P20263)

Gene Name Pou5F1**About Gene Symbol**

Western blot analysis of POU5F1/Oct-4 in F9 lysates using POU5F1/Oct-4 antibody.



Western blot analysis of Oct-4 (3F5) in 20ug and 40ug mESC lysates using POU5F1/Oct-4 antibody.