

Cytokeratin 7 (4G9) Mouse mAb


ABZ4151

Product Overview


Name	Cytokeratin 7 (4G9) Mouse mAb
Accession(Primary)	P08729
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Species Reactivity	Human
Calculated MW	Calculated MW: 51 kDa; Observed MW: 55 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200

Target information (P08729)


Gene Name	KRT7
About Gene Symbol	KRT7 (Keratin 7) is a keratin cytoskeletal protein involved in cytoskeletal organization. It is an important research target in cell biology, cancer metastasis, and cytoskeletal research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Keratin 7(Cterminus) in T47D, HeLa and HepG2 lysates using Keratin 7(Cterminus) antibody.



Immunocytochemistry analysis of Cytokeratin 7 in HeLa cells using Cytokeratin 7 antibody.



Western blot analysis of Cytokeratin 7 (4G9) in HeLa lysates which is positive expression and MCF-7 lysates which is negative expression using Keratin 7(Cterminus) antibody.

