

Hsp70 1A(2D4) Mouse mAb

ABZ4152


Product Overview

Name	Hsp70 1A(2D4) Mouse mAb
Accession(Primary)	P0DMV8
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human
Calculated MW	Calculated MW: 70 kDa; Observed MW: 70 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200


Target information (P0DMV8)

Gene Name HSPA1A

About Gene Symbol



Western blot analysis of Hsp70(Nterminus) in HeLa and HepG2 lysates using Hsp70(Nterminus) antibody.



Immunocytochemistry analysis of Hsp70 1A (2D4)
in HeLa using Hsp70 (Nterminus) antibody.