


Hsp70 1A(2A11) Mouse mAb

ABZ4153


Product Overview

Name	Hsp70 1A(2A11) Mouse mAb
Accession(Primary)	P0DMV8
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Species Reactivity	Human
Calculated MW	Calculated MW: 70 kDa; Observed MW: 70 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20

Target information (P0DMV8)

Gene Name HSPA1A**About Gene Symbol**

Western blot analysis of Hsp70(Cterminus) in HeLa lysates using Hsp70(Cterminus) antibody.



Immunoprecipitation analysis of Hsp70 1A (2A11)
in HeLa lysates using Hsp70 (Cterminus) antibody.