

Integrin beta 3 Mouse mAb


ABZ4154

Product Overview

Name	Integrin beta 3 Mouse mAb
Accession(Primary)	P05106
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human
Calculated MW	Calculated MW: 87 kDa; Observed MW: 110 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P05106)

Gene Name	ITGB3
About Gene Symbol	ITGB3 (Integrin Subunit Beta 3) is an integrin cell adhesion receptor involved in adaptive immunity. It is an important research target in immunology, flow cytometry, and cancer immunotherapy and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Integrin beta 3(Nterminus) in 15ug and 25ug Leukocyte whole lysates using Integrin beta 3(Nterminus) antibody.