

Aurora A (2E3) Mouse mAb

ABZ4155


Product Overview

Name	Aurora A (2E3) Mouse mAb
Accession(Primary)	O14965
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Cat
Calculated MW	Calculated MW: 46 kDa; Observed MW: 46 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200


Target information (O14965)

Gene Name AURKA


About Gene Symbol



Immunofluorescence analysis of Aurora A (2E3) in HeLa cells using Aurora Kinase A antibody(green) and Histone H3.1 (Phospho-Ser1) antibody (red),and DAPI (blue)



Western blot analysis of Aurora Kinase A in SW480 and COS7 lysates using Aurora Kinase A antibody.



Immunocytochemistry analysis of Aurora A (2E3) in HeLa using Aurora Kinase A antibody.